









Rotors

Rotor Selection Guide

	Working Volume (ml)	Number of Tubes	Description	Models SC110/SC110A SPD-Series ISS 110, AES1010, SPD1010	Models SC210A, AES2010, SPD2010	Rotor Cat. No.
MICROCENTRIFUGES TUBES 	1.2 - 1.6	40	1.5 2.0 ml	■		36-072-7
	1.2 - 1.6	64	1.5 - 2.0 ml	■		36-073-8
	1.2 - 1.6	120	1.5 - 2.0 ml	■		36-780-4
	1.2 - 1.6	200	1.5 - 2.0 ml		■	36-075-1
GLASS AND PLASTIC TUBES 	0.3	100	0.5 ml (8 X 29 mm)	■		36-792-155
	0.3	40	0.4 ml (6 x 50 mm)	■		36-792-156
	0.3	100	0.4 ml (6 x 50 mm)	■		36-792-157
	4	20	12 x 75 mm	■		36-073-2
	4	40	1.5 - 2.0 (12 x 75 mm)	■		36-073-3
	4	72	12 x 75 mm	■		36-780-2
	4	200	12 x 75 mm		■	36-075-1*
	8	10	13 x 100 mm	■		36-073-2
	8	32	13 x 100 mm	■		36-781-9
	8	118	13 x 100 mm		■	36-075-1*
	12	48	16 x 125 mm		■	36-075-3*
	10	60	17 x 100 mm, 16 x 100 mm		■	36-075-2*
	15	8	17.5 x 102 mm, 15 ml Corex, tubes	■		36-792-158
	10	8	18 x 100 mm, 17 x 95, 16 x 100	■		36-073-5
	15	48	18 x 125 mm		■	36-075-3*
	24	4	18 x 150 mm	■		36-792-159
24	32	18 x 150 mm, 16 x 125 mm		■	36-792-160*	
30	6	25 x 105 mm, 30 ml Corex, tubes	■		36-792-161	
40	12	28 x 150 mm		■	36-075-7	
CENTRIFUGE TUBES 	12	10	15 ml conical (16 x 120 mm)	■		36-781-5
	12	30	15 ml conical (17 x 120 mm)		■	36-075-2*
	12	52	15 ml conical (16 x 120 mm)		■	36-781-8*
	40	6	50 ml conical (28 x 115 mm)	■		36-784-2
	40	12	50 ml conical (28 x 115 mm)		■	36-075-7
	40	26	50 ml conical (28 x 115 mm)		■	36-781-7*
40	48	50 ml conical (28 x 115 mm)		■	36-792-162*	
FLASKS 	35	8	50 pear-shaped flask	■		36-778-6
	80	4	100 ml pear-shaped flask	■		36-784-3
	80	8	100 ml pear-shaped flask		■	36-792-140
	200	6	250 ml pear-shaped flask		■	36-778-4
	400	4	500 ml pear-shaped flask		■	36-792-163*
	80	8	100 ml recovery flask		■	36-792-164
VIALS 	2	60	12 x 32 mm	■		36-073-1
	3.0	192	1 dram vial (4 ml)		■	36-792-165*
	2.0	60	12 x 40 mm vials	■		36-073-1
	2.4	12	20 x 47 mm v-vials	■		36-792-168
	3	24	1 dram vials (15 x 45 mm), 4 ml	■		36-792-167
	4	12	20 x 60 mm v-vials	■		36-792-168
	5.6	24	18 x 52 mm mini-scintillation vials	■		36-792-169
	16	12	28 x 60 mm scintillation vials	■		36-073-7
16	50	28 x 60 mm scintillation vials		■	36-792-171*	
BOTTLES	250	8	250 ml centrifuge bottles		■	36-792-140
MICROPLATE 	0.3	2	Multiwell Plates (shallow)	■		36-792-25
	0.3	12	Multiwell Plates (shallow)		■	36-789-8
	1 - 2	8	Multiwell Plates (Beckman deepwell)		■	36-789-9
+ CUSTOM ROTORS AVAILABLE						
<i>For SPD-Series concentrators ONLY</i>						
GLASS AND PLASTIC TUBES 	24	6	18 x 150 mm tubes			36-792-92
MICROPLATE 	0.3	6	Multiwell Plates (shallow)			36-792-172
	1 - 2	2	Multiwell Plates (deepwell)			36-792-141

* For use in SC250DDA SpeedVacs with RAD200 adapter.