

A&D WEIGHING

EJ NEWTON Series Compact Balances



1mg Readability



10mg Readability



100mg Readability



170 Shields Court, Unit 2
Matham, ON L3R 9T6
TEL: (905) 476-6880 ext. 236
FAX: (905) 476-1231

A&D
A&D Weighing
...Clearly a Better Value

Get it from your authorized distributor: www.labequip.com



ask@labequip.com

COMPACT BALANCES

The EJ Newton Series

A new model with 1 mg resolution joins the EJ family of compact balances. Affordable and built-to-last, EJ Series balances are used in many applications including:

- Industrial
- Pharmacy
- Education

With thousands of balances in laboratories and production facilities around the world, the EJ Series has provided years of reliable weighing performance and are backed by A&D's industry-leading 5 year warranty.



Affordable Precision

Now students, researchers and industrial users have a convenient and affordable way to make milligram measurements.

The EJ-123 has standard functions usually found in balances costing twice as much:

- Gravity correction
- Accumulation function
- GLP compliant output



EJ-123
(Draft shield standard)

The EJ-123 makes milligram measurements as easy as, well, 1-2-3:

1. Plug-in and warm-up the balance.
2. Level the balance by adjusting the feet.
3. Start weighing.

Options

EJ-02 & EJ-02C



Quick USB interface. Option EJ-02C includes USB cable.

EJ-03



RS-232C (serial) interface.

EJ-07 & EJ-08



Underhook. EJ-07 for EJ-3000 to 6100. EJ-08 for EJ-1500, EJ-2000.

EJ NEWTON SERIES



EJ-11



Breeze break (draft shield) for all models except EJ-123.

EJ-12



Carrying case. All models except EJ-123.

EJ-13



Density determination kit for models EJ-120 to EJ-610.

EJ NEWTON SERIES SPECIFICATIONS

[illegible]

Specifications subject to change without notice.

Units of Measure

[illegible]

Specifications subject to change without notice.

Options

EL-02C	USB interface with cable
EL-03	RS-232C interface
EL-07	Underhook for EL-3000 to EL-6100
EL-08	Underhook for EL-1500, EL-2000
EL-11	Breeze break (Draft Shield)
EL-12	Carrying case (except for EL-123)
EL-13	Density determination kit for EL-120/ EL-200/EL-300/EL-410/EL-610

Accessories

TB:16 ³ C	AC adapter 120V (Replacement)
TB:16 ⁴	AC adapter 240V (Round prongs)



170 Shields Court Unit 2
Markham, ON L3R 9T5
TEL: (905) 475-5880 ext. 226
FAX: (905) 475-1231



LT: EJ 11/ 13

Get it from your authorized distributor: www.labequip.com



ask@labequip.com