

Flash Point by Tag Closed Tester

test method

For flash point determinations of liquids with a viscosity of below 5.5 centistokes (cSt) at 104°F (40°C) or below 9.5cSt at 77°F (25°C), and a flash point below 200°F (93°C) except cut-back asphalts, those liquids which tend to form a surface film under test conditions and materials which contain suspended solids.

tag closed cup flash tester

- Conforms to ASTM D56 and related specifications
- Gas or Electric Heating

Determines flash points of liquid products by the Tag Closed Cup method. Features stepless variable heat control with reference dial for accurate repeat setting of temperature rate of rise per specifications. Also available with gas burner instead of electric heater. Precision machined cover mechanism simultaneously opens slide shutter and applies test flame to sample at the turn of a knob. Includes liquid bath with constant level overflow, brass test cup, plated brass thermometer ferrules and test flame reference bead. Bath and cover mechanism are constructed of plated brass. Heater is enclosed in a cast aluminum base assembly.

ordering information

catalog no.	description
K14601	Tag Closed Cup Flash Tester, 115V 60Hz
K14671	Tag Closed Cup Flash Tester, 220-240V 50/60Hz

accessories

250-000-09F	ASTM 9F Thermometer Range: 20 to 230°F
250-000-09C	ASTM 9C Thermometer Range: -5 to +110°C
250-000-57F	ASTM 57F Thermometer Range: -4 to 122°F
250-000-57C	ASTM 57C Thermometer Range: -20 to +50°C
K14510	Cover Assembly
K14520	Brass Test Cup



K14601 Tag Closed Cup Flash Tester

specifications

Conforms to the Specifications of:
ASTM D56; IP 304; FTM 791-1101

Electrical Requirements

115V 60Hz 1.3A
220-240VAC, 50/60Hz 0.6A

Included Accessories

Brass Test Cup
Cover Assembly

Dimensions wxdxh in.(cm)

10 x 17.5 x 11.5 (25.4 x 44.5 x 29.2)
Net Weight: 20 lbs (9.1kg)

Shipping Information

Shipping Weight: 30 lbs (13.6 kg)
Shipping Dimensions: 2 Cu. ft.



85 Corporate Drive, Holtsville, New York 11742
1-800-878-9070 (in u.s. only) TEL: +1 631 589 3800 FAX: +1 631 589 3815
Email: sales@koehlerinstrument.com www.koehlerinstrument.com

©2021 Koehler Instrument Company, Inc.

Get it from your authorized distributor: www.labequip.com



ask@labequip.com