

T-Edge 11

The only gravity autoclave upgradable to pre/post vacuum by software



Innovative Technology and Design meets
Maximum Capacity, to Provide
Enhanced Patient and Operator Safety


Tuttnauer
Innovation · Legacy · Partnership

Get it from your authorized distributor: www.labequip.com

Labequip

ask@labequip.com

Optimal Chamber

- 5 extra-long trays for wrapped / non-wrapped loads
- 3 full dental cassettes + 3 half dental cassettes
- More than 80% chamber capacity
- Convenient and easy loading with door opening to 110°

Touchpad

- User friendly touchpad and control panel. Sterilization at the click of a button

Longer Life

- Only Clean Water is used due to dual water tank system



Monitoring

- Simple monitoring using your smartphone



Closed Door Drying

- keeps door locked during drying which produces bone dry loads.

Capabilities

- Regulatory & Process Capabilities: Monitoring & Traceability
- Internal memory for storing up to 1,000 sterilization cycles
- USB outlet for transferring data to your PC or network
- Ethernet connection for transferring data to your PC or network

Automatic

- Water Refilling and Draining options available



Service



flexibility



Experience




Quality

Get it from your authorized distributor: www.labequip.com



ask@labequip.com



Eco Friendly

- Reduced water consumption by about 40%
- Decreased energy consumption thanks to reduced cycle time
- Lower electricity costs with time-delay settings for machine to work off hours



Enhanced Serviceability

- Easy to service
- Simple to maintain
- Extended product life

T-EDGE11

Chamber Size	18" depth x 11" diameter
Chamber Volume	7 gallons/27 liters
Dimensions (w x h x d)	21.1" x 19.68" x 19.68"
Weight	132 lbs.
Power Supply	120V, Frq 50-60Hz
Number of Trays	5
Connections	1 Ethernet port, 4 USB ports, WIFI, integrated automatic water filling and drain connection
Items Supplied	Drain tube, tray holder

Icon	Cycle Program	Load Type	*T-Edge11 Cycle Times Including Drying (Minutes)
	Unwrapped Instruments 270° F	Unwrapped Instruments	25 Minutes
	Wrapped Pouches 270° F	Wrapped Pouches	55 Minutes
	Unwrapped Delicate 250° F	Unwrapped Delicate Instruments	47 Minutes
	Handpieces 270° F	Handpieces	60 Minutes

Icon	Cycle Program	Load Type	*T-Edge11 Cycle Times Including Drying (Minutes)
	Unwrapped Instruments 273° F	Unwrapped Instruments	26 Minutes
	Wrapped Pouches 273° F	Wrapped Pouches	56 Minutes
	Unwrapped Delicate 250° F	Unwrapped Delicate Instruments	35 Minutes
	Wrapped Delicate 250° F	Wrapped Delicate Instruments	79 Minutes

So much more than just an Autoclave...

Your skills are needed with your patients and your time is of the essence... T-Edge11 will free you up.


- ◆ Compliance with strict international standards
- ◆ Clinically proven sterilization capabilities
- ◆ High-resolution tactile touchpad
- ◆ Innovative original technology
- ◆ Simple-to-use interface
- ◆ Chic design
- ◆ Great user experience

Maximum Sterilization. Minimum Costs.

Advanced Barcode Traceability



Receive barcodes with sterilization data (including expiration date).



Apply barcode label to each pouch.



Enhance your patients' safety

- Multiple traceability options
- Pre-programmed test cycles
- Protect your clinic's reputation
- Compliance with strict safety standards

For Your T-Edge Chamber Clean:

CLEAN & SHINE Your Chamber

Tablets CB0030 For T-Edge Only




T-Clean Chamber Shine SH0006 500 ml with sprayer





Tuttnauer USA Co. Ltd.

25 Power Drive, Hauppauge, NY 11788

Tel: +800 624 5836, +631 737 4850,

Fax: +631 737 0720

E-mail: info@tuttnauerUSA.com

www.tuttnauerUSA.com



Tuttnauer

Innovation · Legacy · Partnership

Learn from our experts.

Visit our blog: tuttnauerusa.com/blog