

A&D WEIGHING

EJ NEWTON Series Compact Balances



1mg Readability



10mg Readability



100mg Readability

Labquip

170 Shields Court Unit 2
Warham, ON L3R 9T5
TEL: (905) 476-8880 ext. 226
FAX: (905) 476-1291

A&D

A&D Weighing

...Clearly a Better Value

COMPACT BALANCES

The EJ Newton Series

A new model with 1 mg resolution joins the EJ family of compact balances. Affordable and built-to-last, EJ Series balances are used in many applications including:

Industrial
Pharmacy
Education

With thousands of balances in laboratories and production facilities around the world, the EJ Series has provided years of reliable weighing performance and are backed by A&D's industry-leading 5 year warranty.

□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □
	□ □ □ □ □ □ □ □ □ □	□ □ □ □ □ □ □ □ □ □

Affordable Precision

Now students, researchers and industrial users have a convenient and affordable way to make milligram measurements.

The EJ-123 has standard functions usually found in balances costing twice as much:

- Gravity correction
- Accumulation function
- GLP compliant output



The EJ-123 makes milligram measurements as easy as, well, 1-2-3:

1. Plug-in and warm-up the balance.
2. Level the balance by adjusting the feet.
3. Start weighing.

Options

EJ-02 & EJ-02C



Quick USB interface. Option EJ-02C includes USB cable.

EJ-03



RS-232C (serial) interface.


EJ-07 & EJ-08



Underhook. EJ-07 for EJ-3000 to 6100. EJ-08 for EJ-1500, EJ-2000.

EJ NEWTON SERIES

127 x 140 mm stainless steel weighing pan for EJ-1500 to EJ-6100.
 110 mm round pan for EJ-610 & below.



EJ-11



Breeze break (draft shield) for all models except EJ-123.

EJ-12



Carrying case. All models except EJ-123.

EJ-13



Density determination kit for models EJ-120 to EJ-610.

