

Mini Vortex Mixers



Four Mini Vortex Mixer models for gentle to high-speed mixing are available. Analog model offers variable speed control, and digital model allows input of the exact speed and time to achieve reproducible results. Pulsing model features a unique pulse action that reduces heat generation, while providing more effective mixing and cell disruption. Fixed speed model offers one-touch, high-speed mixing.

- Select from Two Modes of Operation—Touch or Continuous
- Variable Speed Models Offer Low RPM Start Up to High-Speed Mixing
- Fixed Speed Models Provide Full RPM for Vigorous Mixing



Mini Vortex Mixers

- Sturdy design
- Fixed high speed mixing
- Touch mode operation

Built sturdy to provide stable and reliable vortexing action. Starts mixing when the cup head is pressed down. Speed is fixed at full rpm to provide vigorous vortexing of samples.

Operating Conditions:

Unit can be run in conditions from 4 to 40°C (39 to 104°F), 20% to 85% relative humidity, non-condensing.

Ordering Information:

Units include a 3-wire cord and plug. Includes both cup head and 3" head with cover.

This product includes:

Description	Item Number
Cup Head	30400235
3" (7.6 cm) Head	30400236
3" (7.6 cm) Rubber Head Cover	30400237



Description	Model	Item Number
Fixed Speed Vortex Mixer	VXMNFS	30392108

Specifications	
Speed Range 120V	3200 rpm
Orbit	4.9 mm (0.19")
Controls	None
Duty Rating	Intermittent duty
Dimensions (L x W x H)	8.3" x 4.8" x 6.5" (21.1 x 12.2 x 16.5 cm)
Ship Weight	10 lbs (4.5 kg)



- Sturdy design
- Variable, analog speed control
- Continuous or touch mode operation

Built sturdy to provide stable and reliable vortexing action. Control allows low rpm startup for gentle shaking or high speed mixing for vigorous vortexing of samples. Two modes of operation; continuous mode when using accessory attachments or touch mode which activates mixing when depressing the cup head.

Operating Conditions:

Unit can be run in conditions from 4 to 40°C (39 to 104°F), 20% to 85% relative humidity, non-condensing.

Ordering Information:

Units include a 3-wire cord and plug. Includes both cup head and 3" head with cover. Additional accessories can be found on pages 52-53.

This product includes:

Description	Item Number
Cup Head	30400235
3" (7.6 cm) Head	30400236
3" (7.6 cm) Rubber Head Cover	30400237



Description	Model	Item Number
Analog Vortex Mixer	VXMNAL	30392115

Specifications	
Speed Range* 120V	300 to 3200 rpm
Orbit	4.9 mm (0.19")
Controls	Auto/Off/On Rocker Switch, Speed Knob: Variable 1 to 10 Dial Marks
Duty Rating	Intermittent duty
Dimensions (L x W x H)	8.3" x 4.8" x 6.5" (21.1 x 12.2 x 16.5 cm)
Ship Weight	10 lbs (4.5 kg)



* Maximum speed will vary depending on accessory used.

Mini Vortex Mixer Accessories

Micro-Tube Holder

Mixes (48) 0.25 to 2 mL micro-tubes.
Requires Insert Retainer.

Description	Item Number
Micro-Tube Holder (2 pack)	30400232
Insert Retainer	30400227



9 to 13 mm Tube Holder

Ideal for mixing 5 mL culture tubes and micro-vials.
Requires Insert Retainer.

Description	Item Number
9 to 13 mm Tube Holder (2 pack)	30400229
Insert Retainer	30400227



14 to 19 mm Tube Holder

Ideal for mixing up to (8) 15 mL centrifuge tubes.
Requires Insert Retainer.

Description	Item Number
14 to 19 mm Tube Holder (2 pack)	30400230
Insert Retainer	30400227



20 to 25 mm Tube Holder

Ideal for mixing up to (8) 50 mL centrifuge tubes.
Requires Insert Retainer

Description	Item Number
20 to 25 mm Tube Holder (2 pack)	30400231
Insert Retainer	30400227



Vessel Harness

Mixes Erlenmeyer flasks and media bottles.
Requires Insert Retainer.

Description	Item Number
Vessel Harness (2 pack)	30400228
Insert Retainer	30400227



Cup Head

Designed for mixing 1 tube at a time.

Description	Item Number
Cup Head	30400235



3" Rubber Head Cover and 3" Head

Designed for mixing irregular shaped objects.

Description	Item Number
3" (7.6 cm) Rubber Head Cover	30400196
3" (7.6 cm) Head	30400236
3" (7.6 cm) Head with Rubber Cover	30400237



Microplate Holder

Ideal for mixing 96-well plates or deep well blocks.
Requires Insert Retainer.

Description	Item Number
Microplate Holder (2 pack)	30400233
Insert Retainer	30400227



Flat Foam Insert

Ideal for custom applications. Can be cut or drilled to fit your specifications. Requires Insert Retainer.

Description	Item Number
Flat Foam Insert (2 pack)	30400234
Insert Retainer	30400227



Microwell Kit with Retainer

Includes:
• 2 Microplate Holders
• 1 Insert Retainer

Description	Item Number
Microwell Kit with Retainer	30400199



Microtube Kit with /Retainer

Includes:
• 2 Micro-Tube Holders
• 1 Insert Retainer

Description	Item Number
Microtube Kit with /Retainer	30400200



Foam Insert Variety Pack

Includes:
• (1) 9 to 13 mm Tube Holder
• (1) 14 to 19 mm Tube Holder
• (1) 20 to 25 mm Tube Holder
• 1 Flat Foam Insert

Description	Item Number
Foam Insert Variety Pack	30400201



Micro-Tube and Microplate Holder Kit

Includes:
• 1 Micro-Tube Holder
• 1 Microplate Holder
• 1 Insert Retainer

Description	Item Number
Micro-Tube & Microplate Holder Kit	30400198



Tube Holder Kit

Includes:
• (1) 9 to 13 mm Tube Holder
• (1) 14 to 19 mm Tube Holder
• (1) 20 to 25 mm Tube Holder
• 1 Flat Foam Insert
• 2 Vessel Harnesses
• 1 Insert Retainer

Description	Item Number
Tube Holder Kit	30400197



OHAUS Vortex Mixer accessories are interchangeable with other manufacturers' models. Call 1-800-672-7722 for assistance.

Accessories

Mini Vortex Mixer Accessories

Single Tube Holder

Single tube, hands free mixing. Easily attaches to the top of any Vortex Mixer with the use of a strong magnetic base. Accepts tubes from 2.5 to 4.5" (6.4 to 11.4 cm) in length. Minimum tube diameter of 0.75" (19 mm).

Description	Item Number
Single Tube Holder	30400202



Adapter for Vortex-Genie® Mixer*

Adapter plate easily adheres to the Vortex-Genie® Mixer housing so Single Tube Holder (sold separately) can be attached.

Description	Item Number
Adapter plate	30400203



* The Vortex-Genie® Mixer is a registered trademark of Scientific Industries, Inc.

Single Tube Holder Accessories

(Requires Single Tube Holder)

0.5 mL Micro-Tube Holder

Mixes (24) 0.5 mL micro-tubes.
For use with Single Tube Holder.

Description	Item Number
Micro-Tube Holder (0.5 mL)	30400204



1.5 mL to 2.0 mL Micro-Tube Holder

Mixes (18) 1.5 to 2 mL micro-tubes.
For use with Single Tube Holder.

Description	Item Number
Microtube Holder (1.5 mL to 2.0 mL)	30400205



Ampule Tube Holder

Mixes up to 4 storage vials and test tubes.

Description	Item Number
15 to 17 mm Ampule Tube Holder	30400207
10 to 17 mm Ampule Tube Holder	30400208



Stainless Steel Microtube Holder

Mixes up to (12) 1.5 to 5 mL tubes.
Stainless steel construction.

Description	Item Number
Microtube Holder	30400206

