



Desktop High Energy Ball Mill

VQ-N ball mills may be referred to as shaker mills, mixer mills or high-energy ball mills. They can be used for both wet and dry grinding as well as to mix a variety of solids, suspended liquids and pastes. This unit is ideal for preparing small batches of your sample efficiently.



Features

- Quick, efficient pulverization and homogenization
- Short grinding time and high sample throughput
- Large range of materials available for grinding media
- Leak-proof jar for wet grinding

Model	VQ-N
Electrical Requirements	110V-220V 50/60Hz single phase, 180 watts
Jars Capacity	50ml – 80ml
Oscillation Speed	1,200 rpm
Timing Range	1-999 seconds
Discharging Granularity	Down to 0.1 micron
Grinding Jar Included	One 80ml stainless steel jar with stainless steel balls and o-ring. Additional grinding media in a variety of materials is sold separately.