

MELTING POINT OF PETROLEUM WAX (COOLING CURVE)



K17500 Wax Melting Point Apparatus

Test Method

Periodic temperature measurements are taken of a sample of molten wax as it is cooled in an air bath. When the wax solidifies, a plateau in the cooling curve occurs, indicating the melting point (cooling curve) of the sample.

Wax Melting Point Apparatus

- Conforms to ASTM D87 and related specifications

Cools molten wax samples in accordance with ASTM and related specifications. Consists of nickel-plated air and water bath assembly with removable cover. Supports test tube in a vertical position in the air bath.

Specifications

Conforms to the specifications of:

ASTM D87; IP 55; ISO 3841; DIN 51570; FTM 791-1402; NF T 60-114

Included Accessories

Test Tube, Thermometer Holders (2)

Dimensions dia.xh,in.(cm) 5½x7 (14x18)

Net Weight 4 lbs (1.8kg)

Shipping Information

Shipping Weight: 6 lbs (2.7kg)

Dimensions: 0.7 Cu. ft.

Ordering Information

Catalog No.		Order Qty
K17500	Wax Melting Point Apparatus	1
Accessories		
250-000-14F	ASTM 14F Thermometer Range: 100 to 180°F	2
250-000-14C	ASTM 14C Thermometer Range: 38 to 82°C	
K175-0-8	Test Tube, 25x100mm	

For NIST traceable certified thermometers, please refer to the ASTM Thermometer section on pages 184 through 191.