



An anti vibration table is designed for cooperation with analytical and laboratory balance or control scales. Its construction provides stable operation during precise weighing processes.

Working table has a separate basis and a top made in painted MDF technology. The top has an opening for stone plate which is a stabilizing element of the anti vibration table. Table construction is made of mild steel profiles with adjustable feet for height regulation. A stone plate is supported by an independent basis for improvement of stability.



Anti vibration table components as well as the top can be manufactured in stainless steel technology (acid-proof), and table dimensions are adjustable to customer requirements.

Anti vibration table is accessible in three options:

SAL / C – anti vibration table in mild steel technology for laboratory balances

SAL / H – anti vibration table in stainless steel technology for laboratory balances

Dimensions:



Technical data:

Marble plate:		Working table:	
Length	410 mm	Length	770 mm
Width	270 mm	Width	650 mm
Height	115 mm	Height	815 mm
Weight	~34 kg	Weight	~21 kg

RADWAG USA L.L.C.

19599 NE 10th Ave., #E, • North Miami Beach, FL, 33179 • USA

Phone: 1-305-651-3522 • Fax: 1-305-651-3523 • <http://radwagusa.com> • e-mail: office@radwagusa.com